

MIMSAL

LUXIFLEX SENSOR PLUS II

Examination and **advanced exploration luminaire.**

Version

Luxiflex sensor Plus II
50.000 lux at 50 cm



High efficiency

The powerful LED light source offers exceptional light quality and a lifespan of more than 50,000 hours.

The Plus II version has a high performance of the color Ra 97% that depicts the exceptional quality of white light for an excellent distinction of fabrics, reproducing all visible colors, naturally and precisely.

Contactless switch

Allows to switch the lamp on and off without contact. Adjustable light intensity (100% - 50% - 0%).

Flexibility

Flexible, height-adjustable and highly precise arm that allows the luminaire to be adjusted to any position.

Design & handling

Designed with clean and soft shapes which facilitate hygiene and cleanliness in medical environments.

Its compact head and the switch located at the base of the lamp make this model an ideal luminaire for the most demanding examinations procedures.

Mounting options

The standard version includes AH table clamp.

As for accessories, we offer a wide selection of supports such as wall supports, wall extension arms, trolley stands and rail brackets.



LUXIFLEX SENSOR PLUS II

Technical Details	(REF. LX2324SP)
Light source	LED 12W
Number of LEDs	4
Useful life	50.000 hours
Illumination at 50 cm	50.000 lux
Color temperature	4.500 °K
Color rendering	>97%
Total length	935 mm
Color	White
Dimmable	YES

ACCESSORIES



Wall B Support
REF. B0101800



Ah Table Clamp
REF. B010190



Rail Plus Bracket
REF. 82026



Extension Arm
REF. B0101S0



Trolley stand
4 kg REF. 09655
8,8 kg REF. 09651