

## MEDICAL FURNITURE



We develop and produce specialized medical furniture for private medical practices. Moreover, we offer you solutions in organizing your private medical practice in order to transform it into a pleasant working place for you and your patients.

Our goal is to create furniture that fits to the chosen environment ergonomically as well as esthetically.



## Examination table with stainless steel frame

Order number: PC-190 TPI-1

### Technical description:

- Made in Romania and registered as 1st class device by the Health Ministry
- Material: square stainless steel pipe (type AIS304 / 304 L), sterilizable with disinfection solutions.
- Upholstery color: various colors available
- Mattress: 6 cm polyurethane sponge upholstered with anti-bacterial eco-leather
- Adjustable headrest in 10 positions with 0-45 degrees inclination, with 120 kg loading capacity.
- Floor protection legs, fitted with non-slip puffs
- Size (LxWxH): 190x65x80 cm
- Stainless steel stabilizers by 2 screws type DIN 967 at the 4 support legs (l = 175.5cm)
- Total thickness of mattress: 8 cm
- Netto weight: 32.5 kg
- Maximum load: 250 kg



### Standard colors



B365\_1044\_AS



B365\_6263\_AS



B365\_9001\_AS

## Optional accessories:

### SS-P-0 - Infusion stand

#### Technical description:

- Material: stainless steel
- Adjustable in height: 100-200 cm
- 4 hooks for holding infusion bags (Max load: approx. 6 kg / hook)
- Netto weight: approx. 0.4kg



### SM-P-0-Disposable glove holder

#### Technical description:

- Material: stainless steel
- Width: 12.5 cm
- Length: 12 cm
- Height: 23 cm
- Netto weight: approx. 0.3 kg



### SP-P-0-Paper roll holder

#### Technical description:

- Material: stainless steel
- Width: 23.5 cm
- Length: 63 cm
- Height: 4 cm
- Netto weight: approx. 0.45 kg



### SSP-0 - Goepel leg support

#### Technical description:

- Material: stainless steel
- Adjustable height
- Upholstered with eco-leather



### SC-0 Steps

#### Technical description:

- Material: rectangular stainless steel frame and seethrough 10 cm plexiglass
- Size (LxWxH): 46x60x34 cm.
- 2 steps.



## Examination table with built in furniture

Order number: **PC-190 MP-0**

### Technical description:

- Material: square stainless steel pipe (type AIS304 / 304 L), sterilizable with disinfecting solutions.
- Upholstery color: various colors available
- Mattress: 6 cm polyurethane sponge upholstered with anti-bacterial eco-leather
- Adjustable headrest in 10 positions with 0-45 degrees inclination, with 120 kg loading capacity.
- Floor protection legs, fitted with non-slip puffs
- Furniture made of 18 mm thick white melamine fiberboard with 4 sections all on wheels. Wheel diam. 50 mm. 3 sections with shelves and 1 with drawers. Can be used separately from the examination table. Size: 190x65x78 cm
- Stainless steel stabilizers by 2 screws type DIN 967 at the 4 support legs (l = 175.5 cm)
- Total thickness of mattress: 8 cm
- Netto weight: 100 kg
- Maximum load: 250 kg



### Optional accessories:

- SS-P-0 Infusion stand
- SP-P-0 Paper roll holder
- SM-P-0 Disposable glove holder
- SSP-0 Goepel leg support
- SC-0 Stairs



## Massage bed

Order number: **PC-190 TP-0**

### Technical description:

- Mattress 6 cm polyurethane sponge upholstered with anti-bacterial eco-leather
- Rectangular stainless steel frame
- The head-rest adjustable in 10 positions with 0-45 degrees inclination and Head rest supporting weight: 120kg
- With facial hole
- Maximum load: 250 kg
- Size (LxWxH): 190x65x75 cm
- Netto weight: 30 kg;
- Various colors available

## Examination table with painted steel frame

Order number: **PC-190 TPO-1**

### Technical description:

- Material: electrostatically white painted square steel pipe
- Various colors available
- Mattress: 6cm polyurethane sponge upholstered with anti-bacterial eco-leather
- Adjustable headrest in 10 positions with 0-45 degrees inclination, with 120kg loading capacity.
- Legs with floor protection, fitted with non-slip puffs
- Size: 190x65x80 cm
- Stabilizers made of steel by means of 2 screws type DIN 967 to the 4 support legs (l = 175.5 cm)
- Mattress: 8 cm
- Netto weight: 32.5 kg
- Maximum load: 250 kg



### Optional accessories:

- SS-P-0 Infusion stand
- SP-P-0 Paper roll holder
- SM-P-0 Disposable glove holder
- SSP-0 Goepel leg support
- SC-0 Stairs



## Flexible endoscope support with 2 tubes

**Order number: SE-F-0**

### Technical description:

- Individually removable tubes made of plexiglass material
- Equipped with 1 full end tube for disinfection, 1 free end tube for drying/storing
- Tube support plate size: 240x100x10mm (LxWxH)
- Tube size: 500 mm
- Diam. interior 36 mm
- Netto wieght: 800 g
- Possibility of mounting on the wall or medical trolley
- Optional accessories: camera head support through C mount coupler



## Flexible endoscope support with 3 tubes

**Order number: SE-F-1**

### Technical description:

- Individually removable tubes made of plexiglass material
- Equipped with 2 full end tubes for disinfection/washing, 1 open end for drying/storage
- Tube support plate size: 240x100x10mm (LxWxH)
- Tube size: 500 mm
- Diam. interior 36 mm
- Netto weight: 950 g
- Possibility of mounting on the wall or medical trolley
- Optional accessories: camera head support through C mount coupler



## Flexible endoscope support with 4 tubes

**Order number: SE-R-0**

### Technical Description:

- Individual removable tubes, made of plexiglass material
- Equipped with 2 full end tubes for disinfection, 2 open end tubes for drying/storage
- Tube support plate size: 240x100x10 mm (LxWxH)
- Tube size: 300 mm
- Diam. interior 21 mm
- Netto weight: 650 g
- Possibility of mounting on the wall or medical trolley



## Flexible endoscope support with 5 tubes

**Order number: SE-R-1**

### Technical Description:

- Individual removable tubes, made of plexiglass material
- Equipped with 2 full end tubes for disinfection, 3 open end tubes for drying/storage
- Tube support plate size: 240x100x10 mm (LxWxH)
- Tube size: 300 mm
- Diam. interior: 25 mm
- Netto weight: 750 g
- Possibility of mounting on the wall or medical trolley



## 3-section medical screen

**Order number: P-170 I2-0**

Technical description:

- Material: stainless steel, non-woven polypropylene canvas
- Foldable
- Equipped with 2 wheels
- Width: 60 cm (folded) / 190 cm (open)
- Height: 175 cm
- Netto weight: 10.5 kg



## 4-section medical screen

**Order number: P-180 I8-0**

Technical description:

- Material: stainless steel, non-woven polypropylene canvas
- Equipped with 8 wheels, 2 wheels on each section
- Width: 50 cm (folded) / 195 cm (open)
- Height: 181 cm
- Netto weight: 10.5 kg



## Negatoscope

**Order number: Neg-0**

Technical description:

- Material: aluminum frame
- Lighting with LED bulbs
- Dimmable light intensity controlled by remote control
- Inox clip for film fixing
- Horizontal mounting
- Size (LxWxH): 90x60x48 cm
- Screen size: 74x40 cm
- Netto weight: 7.6 kg
- Power supply: 220 V
- Dimmable light:
  - 25%: 1500 lux
  - 50%: 3000 lux
  - 75%: 4500 lux
  - 100%: 6000 lux



## Changing table with shelf

**Order number: Minf-P-0**

### Technical description:

- Material: 18 mm melamine fiberboard
- With sliding side shelf on the rail for towels and diapers
- Drawer
- 4 wheels (2 lockable)
- Size (LxWxH): 54x70x93 cm - without wheels
- Netto weight 25 kg
- Drawer: 53 / 50 cm
- Sliding shelf: 55 / 65 cm



## Hospital bassinet

**Order number: PNN-O**

### Technical description:

- Materials: Stainless steel frame, acrylic, seethrough crib
- Adjustable height: 100-115 cm
- Base with 4 wheels (2 lockable);
- Mattress: 5 cm thick sponge upholstered with eco leather
- Size (LxW): 80x100 cm;
- Crib size (LxWxH): 80x45x28 cm
- Netto weight: 4.5 kg;
- Trendelenburg mechanism - inclination  $\pm 15^\circ$

## Breastfeeding chair

**Order number: SC-AI-0**

### Technical specifications:

- Size (LxWxH): 60x55x93 cm
- Netto weight: 12 kg
- Maximum load: 150 kg
- Various colors available
- Antibacterial eco-leather made of PVC, phthalan free, non-flammable
- Seat with 180 mm whole to ensure comfort during breastfeeding



## Medical work table

**Order number: MIRB-0**

### Technical description:

- Material: stainless steel
- Equipped with 4 rotating wheels, of which 2 lockable
- Side tray for storing containers
- Waste bowl holder, interior diam. 23 cm
- Size (LxWxH): 72.5x40x90 cm
- Distance between the two shelves: 58 cm
- Netto weight: 12 kg
- Maximum load: 20 kg



## Bandage table

**Order number: MP-1**

### Technical description:

- Rectangular pipe frame 25x25x1,5 mm, electrostatically painted.
- Drawer material: melamine fiberboard.
- Shelf material: 1,5mm plywood
- Size (LxWxH): 44x80x84 cm.
- Width with handle is 87 cm.
- Netto weight: 16 kg
- Includes wastebowl holder diam.: 75 mm
- Equipped with 4 mobile wheels and 360 degree rotation





## Radio Surg Cart


**Order number: CA-R-0**

### Technical description:

- Equipped with 4 wheels of which 2 lockable
- Welded base made of aluminum pipe, painted in black
- Plexiglass shelves: 39x42 cm
- 18 mm thick melamine fiberboard drawer with silent rail
- Stainless steel basket for tools and accessories
- Size (LxWxH): 50x50x110 cm
- Netto weight: 20 kg

### Available colors

snow		WH10
plum		4H01
pine		6C01
terra		8C01



## Medical endoscopy cart

**Order number: CM-A2-0**


### Technical description:

- Equipped with 4 wheels out of which 2 are lockable
- Made of aluminum
- 2 height adjustable plexiglass shelves
- Key lock drawer
- Monitor holder with VESA monitor clamp
- 4 power outlets with power switch
- Protective spiral tube for power cables
- Size (LxWxH): 50x50x160 cm
- Netto weight: 14 kg
- Various frame color
- Various shelf color

### Optional accessories:

- Flexible endoscope holder with 2 tubes
- Flexible endoscope holder with 3 pipes and optional camera holder
- Rigid endoscope holder with 5 tubes
- Rigid endoscope support with 4 tubes

### Available colors

snow		WH10
plum		4H01
pine		6C01
terra		8C01



## Medical cart

**Order number: CM-I2-0**





### Technical description:

- Equipped with 4 wheels out of which 2 are lockable
- Made of stainless steel
- 2 height adjustable shelf with a click system
- Key lock drawer
- Monitor holder with VESA monitor clamp
- 4 socket connector and switch for powering medical devices
- Protective spiral tube for power cables
- Size (LxWxH): 50x50x160 cm
- Netto weight: 14 kg
- Shelf color: at the beneficiary's request

### Optional accessories:

- Flexible endoscope holder with 2 tubes
- Flexible endoscope holder with 3 pipes and optional camera holder
- Rigid endoscope holder with 5 tubes
- Rigid endoscope support with 4 tubes

### Available colors

snow		WH10
plum		4H01
pine		6C01
terra		8C01





## EROND light UVB 311-108

Order number: UVB-108-P-0

### Technical description:

- Matte stainless steel sheet base 1 mm, shiny stainless steel sheet 0.5 (mirror)
- The protective cast made of 6mm thick plexiglass
- Support made of stainless steel pipe D20x1.25 mm and D16x1.5 mm
- 360 °Rotatable with stabilizing possibility
- Key operated power switch (for authorized personnel only)
- Size (WxHxD): 57x37x5,7 cm (without support); 57x37x21 cm (with support);
- Active surface 400x280 mm
- Weight (with power cord): approx. 5 kg (without support), approx. 6 kg (with support)
- Voltage: 230 V ~ / 50 Hz - 60 Hz
- UV tubes: Philips PL-L 36W / 01 / 4P
- Emission spectrum: 311nm
- Installed electric power: 110W
- Tube electric power: 108W (3 UVB tubes of 36W)
- Adjustable time range: 1 sec. – 15 min.
- Overwriting time resolution: 1 sec.
- Radiation intensity 10 cm from plexiglass: 6 mW / cm<sup>2</sup>
- Radiation intensity 15 cm from plexiglass: 3.8 mW / cm<sup>2</sup>
- Radiation intensity 20 cm from plexiglass: 2.5 mW / cm<sup>2</sup>



## 3 wheel infusion stand with welded base

Order number: SP-200 IP-3-1

### Technical description:

- Material: shiny stainless steel.
- Base diameter: 65 cm.
- Netto weight: approx. 3.7 kg.
- 4 hooks (min. 6 kg / hook).
- Adjustable height: 112-195 cm.
- Welded square pipe base with 3 wheels of which 1 wheel lockable

## Infusion stand with bended base

Order number: SP-200 IP-3-2

### Technical description:

- Material: shiny stainless steel.
- Base diameter: 65 cm.
- Netto weight: approx. 3.7 kg.
- 4 hooks (min. 6 kg / hook).
- Adjustable height: 112-195 cm.
- Bended square pipe base with 3 wheels of which 1 wheel lockable



## Waste bowl with 2 bowls

**Order number: B-24-IP3-6**

### Technical description:

- Material: Shiny, tubular stainless steel frame min. D25 and D20
- 25x25 mm bended stainless steel base, with 3 D75 wheels, of which 1 lockable wheel
- Diameter. base: 65 cm
- Adjustable height: 73-120 cm
- Netto weight: 3.5 kg
- Includes 2 stainless steel bowls, ø33 cm, 6 liters
- Optional: the diameter and volume of the bowls can be chosen because the support is adjustable



## 5-wheel base waste bowl

**Order number: B-24-IP3-6**

### Technical description:

- Material: shiny tubular stainless steel D25 and D20
- Metal base on 5 wheels.
- Base diam.: 65 cm
- Adjustable height: 73-120 cm
- Netto weight: 3 kg
- Includes 1 stainless steel bowl, ø33 cm, 6 liters
- Optional: the diameter and volume of the bowls can be chosen

## Manual pump for liquid Nitrogen

**Order number: PMC-1**

### Technical description:

- Self-priming pump specially designed for liquid nitrogen cylinders, without losses.
- Height: 700 mm.
- Netto weight: 1 kg.
- Suction depth: 500 mm
- Flow rate: 1.5l / min
- Material: shiny stainless steel, tubular.
- The polyamide plug ensures the pressure in the container.
- Bronze filter prevents the penetration of ice crystals into the surgical unit, ensuring a smooth operation



## Liquid Nitrogen scoop

**Order number: SR-TP-0**

### Technical description:

- Material: shiny stainless steel.
- Extraction volume approximately 100 ml of liquid nitrogen.
- Plastic covered handle.
- Extraction depth: 300 mm.
- Total length: 460 mm.
- Netto weight: 90g

## Germicidal lamp 2x55W

**Order number: LG-2x55 M-1**

### Technical description:

- UV-C lamps: 2 55W UV-C tubes
- Size (LxWxH): 24x13x57,5 cm
- Power supply: 230 V ~ / 50 Hz - 60 Hz
- Emission spectrum: 254nm
- Power cable length: 3 m
- Timer: 2 h
- Electrical protection class: I, IP20
- Ozone -free
- During the operation of the device, personnel access is not allowed
- Support of choice



## Germicidal lamp 3x18W

**Order number: LG-165 M-1**

### Technical description:

- UV-C lamps: 3 18W UV-C tubes
- Size (LxWxH): 25,5x12x26 cm
- Power supply: 230 V ~ / 50 Hz - 60 Hz
- Emission spectrum: 254nm
- Power cable length: 3 m
- Timer: 2 h
- Electrical protection class: I, IP20
- Ozone -free
- During the operation of the device, personnel access is not allowed
- Support of choice

## LAMP SUPPORTS

### SL-P-0 Wall hanger

#### Technical description:

- Stainless steel pipe and sheet (DIN 1.4301)
- Mountable on the wall with 4 screws
- Size (LxWxH): 18x10x21 cm
- Netto weight: 0.4 kg



### SM-LG-0 Mobile stand with wheels

#### Technical description:

- Stainless steel round and rectangular pipe (DIN 1.4301)
- Base diameter: 57 cm
- Adjustable height: 138-188 cm.
- Netto weight: 3.5kg
- Bended stainless steel base with 3 wheel out of which 1 lockable
- Wheel diam: 75 mm



### TL-IP3-1 Desktop stand

#### Technical description:

- Stainless steel round and rectangular pipe (DIN 1.4301).
- Base diameter: 350 mm
- Height: 440 mm
- Netto weight: 0.7 kg





RO - 540236, Targu Mures  
Gh. Doja Street 177, Prodcomplex, C1 building  
+4 0724 044 076 | office@erond.ro | www.erond.ro

