



MedGyn

Innovative Products for
Ob/Gyn Professionals

MedGyn (Office Use) Endoscopic (Hysteroscopic) Cannula

MedGyn's Endoscopic (Hysteroscopic) Cannula is a single-use cannula for uterine observation pre- and post-surgery.

Features & Diagnostic Benefits

- Sterile, single-use, highly-durable lightweight, water resistant 7mm cannula
- High resolution eagle eye optical system
- Provides viewing of entire uterine cavity
- Anatomical curve designed cannula with dual function:
Fluid infusion and suction
- Cannula has connections for light source and 8" monitor
- Ideal diagnostic device for: Cases of abnormal uterine bleeding, infertility, uterine polyps and fibroids, suspected uterine or cervical cancer, repeated miscarriages and uterine malformations
- Workstation includes cannula, monitor and foldable stand.
Up-right rolling stand is optional



Monitor with Stand



Monitor without Stand




Endoscopic (Hysteroscopic) Cannula





Image Capture Software with Laptop



Therapeutic For

- Removal of endometrial polyps
- Removal of small fibroids
- Removal of adhesions (Asherman's Syndrome)
- IUD setup and removal
- Dilatation and Curettage

Specifications

Illumination > 500 LX

Field of View > 60 Degree

Resolution > 201 p/mm

Working Distance 170 mm

018145 MedGyn Endoscopic (Hysteroscopic) Cannula

018146 MedGyn Endoscopic (Hysteroscopic) Cannula and Monitor with Stand

018147 MedGyn Endoscopic (Hysteroscopic) Cannula and Monitor without Stand

018148 MedGyn Endoscopic (Hysteroscopic) Cannula Image Capture Software with Laptop

018149 MedGyn Endoscopic (Hysteroscopic) Cannula Image Capture Software

TRUSTED PRODUCTS FOR PREMIER WOMEN'S HEALTHCARE

MEDGYN PRODUCTS, INC.

100 W. Industrial Rd. Addison, IL 60101 USA

t: +1 630.627.4105 toll-free: 800.451.9667

f: +1 630.627.0127

e: info@medgyn.com

w: medgyn.com

ENDOSCOP-PF2

