

**KERNEL** 科诺



# KN-601A

Multi-parameter Patient Monitor





## Functions and Features

- 7"(16:9)narrow screen TFT true color high definition LCD , small volume, light weight
- Real-time display the waveform of ECG, SPO2 and RESP. At the same time, display the digital data of HR/PR, TEMP, RESP, SPO2, NIBP etc. ETCO2 is Optional.
- Three monitoring modes: monitor, operation and diagnosis
- Blood pressure measurement has three modes: manual, cycle and fast measurement
- Have the freezing and storage function
- Display the latest 400 NIBP measured data
- Optional inner recorder, recording many waveforms and events
- Store 96 hours trend chart
- Three levels alarm, free setting function
- Anti-defibrillation and anti-electrosurgical unit, and can be used simultaneously with the pacemaker
- AC-DC. work time of the battery is over 1 hour
- LAN interface, can connect KN-601 Central Monitoring System any time
- Many fixed modes such as bottom bracket, back bracket etc



## **KERNEL MEDICAL EQUIPMENT CO., LTD.**

Add:Kernel Building, Economic Development District,  
Xuzhou City, Jiangsu Province, China  
Tel : 0086-0516-87732209 87732210 Fax : 0086-0516-87732210  
E-mail:admin@kernelmed.com Http://www.kernelmed.com