

## SM-K-A anesthesia trolley

- stainless steel made
- locking upper drawer
- Ø 125 mm castors, 2 equipped with brakes




**SM-K-A**  
506.100  
Dimensions (l/w/h): 1060x655x915 mm



**SM-K-A1**  
506.200  
Dimension (l/w/h): 760x655x915 mm

## SI instrument tableструментальные


- stainless steel made



**SI-60**  
502.015  
Dimensions (l/w/h):  
1100x650x850 mm



**SI-60-1**  
502.025  
Dimensions (l/w/h):  
1100x650x850 mm



**SI-60-2**  
502.005  
Dimensions (l/w/h):  
1200x675x850 mm



**SI-61**  
502.006  
Dimensions (l/w/h):  
1700x650x850 mm



**SI-61-1**  
502.026  
Dimensions (l/w/h):  
1700x650x850 mm



**SI-61-2**  
502.016  
Dimensions (l/w/h):  
1800x675x850 mm