


FlowZig



The wireless UroFlowmetry



**Divisione
Menfis**

Sede Operativa:

Via della Beverara, 46/D - 40131 Bologna Italy.

Tel +39 051 635 10 74 - Fax +39 051 635 10 81


E-mail: info@menfis.it - Web site: www.menfis.it

MEDICA S.p.A.

FlowZig is the Medica - Menfis Division Wireless Uroflowmeter

Wireless Uroflowmetry with
BLUETOOTH® protocol
Low power consumption
Band 2402-2480 MHz
Max. operative distance 50 m (in open space)
Compatible with any Laptop or Desk Top
working with Windows Operating System
Software included

Analysis
Uroflowmetry
Manual and Automatic Analysis
Siroky maximum and average
flow rate nomograms
Liverpool Nomograms



Flowmeter

Load cell flow trasducer at high resolution
and long term accuracy
Flow range 0 to 100 ml/s
Accuracy 1%
End scale: automatic or manual
Sampling rate: 10Hz
A/D converter: 12 bit
Volume range: 2000 ml



* Technical data can be changed without any notice

Rev. 2.0 Mar'16

**Menfis
Division**

MEDICA S.p.A.
via degli artigiani, 7
41036 MEDOLLA MODENA - ITALY

Divisional Office:
Via della Beverara, 46/D - 40131 Bologna Italy.
Tel +39 051 635 10 74 - Fax +39 051 635 10 81
E-mail: info@menfis.it - Web site: www.menfis.it

CE
0476