



RF13A



Portable High-Definition Digital Cystoscope



The New Standard in Flexible Video Cystoscopy and Reprocessing



Compatible with a wide range of surgical procedures

The large $\varnothing 2.2$ mm channel allows you to use a wide range of highly versatile instruments. Electrosurgical treatments such as hemostasis and hot biopsies can also be performed.



High-definition images clearly observe pathological details

To satisfy demand for less invasive and more efficient urologic endoscopes. Advanced video imaging capability ensures outstanding observation capability with clear, vivid reproduction of internal tissues structures.



Model	Outer diameter	Inside diameter	Bending	Working Length	View Angle	Image Definition
RF13A	4.9mm	2.2mm	U210°/D130°	380mm \pm 10%	90°	HD
	5.5mm	2.6mm				