

# RF14A

Flexible Video Ureteroscope



**U**ltra-light handle, long-term grip without getting tired.

**R**eusable, supports liquid and low temperature plasma sterilization, and reduces the cost of use.

**S**lim outer diameter, better water return effect.



### Optical system


- Field of view 110°
- Depth of field 3–50mm
- Direction of view 0°

### Insertion portion

- Insertion tube OD 7.5Fr/8.4Fr/9.2Fr
- Angulation range U 275° / D 275°
- Inside diameter 3.6Fr
- Working length 680mm

8.4Fr

8.4Fr  
No stone left



275°  
Cover all at one glance



16,000  
Ultra-clear imaging



weseemed  
medical technology

**ZHONGSHAN WESEE MEDITECH CO., LTD.**

Address: 6F, Block C, 2Jianshuo Str. HoChung Village, Nanlang, Zhongshan 528451 P. R. China

TEL: (+86)0760-86960005

Email: info@weseemed.com

www.weseemed.com