





weseemed
medical technology

M110A

VIDEO PROCESSOR



Main Feature



Multiple types of endoscopes can be connected, plug and play, convenient and fast.



Dual HDMI ports can connect two medical monitors at the same time, which is convenient for multiple people to watch the teaching.

Image

■ Display resolution	1280×800
■ White Balance	Auto white balance
■ Power Supply	Open-use
■ Basic Function	Taking photo and video
■ Replays Function	Checking recorded photos or videos
■ Storage Mode	SD card
■ Video Output	Provide HDMI and other signal output methods

Power Supply

■ Voltage	110-240V
■ Frequency	50-60Hz



ZHONGSHAN WESEE MEDITECH CO., LTD.

TEL:(+86)0760-86960005

Email:info@weseemed.com

www.weseemed.com